

## MALLEOFIX ACTIVE

**Ankle orthosis, usable on both sides**



### INDICATIONS

- Chronic ligament instability
- Postoperative care
- Acute ankle injuries



### MODE OF OPERATION

Stabilization and functional mobilization of the OSF / USG, in three phases, thanks to modular structure of MalleoFix Active.

- Acute phase - stabilization and immobilization: orthosis with aluminum profiles, supination and pronation
- mobilization - controlled load: orthosis with coil springs, supination and pronation
- Movement - load: orthosis as stabilizing bandage


### FEATURES

C6TEX, the latest generation of elastic textiles.


The honeycombed carbon fiber structure woven into the elastic material increases stability and performance, at the same time it has an antibacterial effect, regulates climate and is antistatic.

- Intuitive handling and functionality through simple and slender construction
- Supination and pronation can be individually positioned
- Customizable aluminum profiles

### ACCESSORIES



Aluminum Extrusions



Helical springs

**Sizes, item number and ankle circumference \***

S	M	L	XL
1000 02030	1000 02031	1000 02032	1000 02033
< 23 cm	23 - 27 cm	27 - 31 cm	>31 cm

**Aid No. 23.02.04.0005**

\* 3 cm above the ankle