



ACTIVE CARPAL

Orthosis for immobilization

INDICATIONS

- Carpal tunnel syndrome
- Wrist arthritis
- Dislocations of the wrist of postoperative
- Immobilization
- Sprains

FEATURES

- Anatomically preformed, cold deformable aluminum rail
- Customizable
- Modern, breathable material combination

Sizes, item number, back width and page (left)

S	M	L
1000 07084	1000 07085	1000 07086
5,5-7 cm	7-8,5 cm	8,5-10 cm

Sizes, item number, back width and page (right)

S	M	L
1000 07080	1000 07081	1000 07082
5,5-7 cm	7-8,5 cm	8,5-10 cm

MODE OF OPERATION

- Immobilization in a functional position
- Fixation of each finger and the thumb
- Joint capsule release
- Protection of the carpal carpal tunnel area

WIDTH BACK OF HAND



Aid No. 23.07.02.3013