

## ERGOLOC

### Symphysis belt / Pelvic ring orthosis



#### INDICATIONS

- ISG syndrom
- Pubic symphysis disruption
- Loosening of symphyses

#### MODE OF OPERATION

- Extension of pelvis
- Stabilization and relief
- Trochanter supports

#### FEATURES

- Completely new in terms of functionality and therapeutic use, by means of patented technology.
- With the Ergoloc, for the first time a pelvic ring orthosis is available. Offers individualized support due to its low height and the proven pulley block principle. Does not restrict patient's movements.
- Through lateral, ergonomically formed pads with variable positioning, the connection between the non-elastic belt material and the patented locking mechanism allows a marked relief and stabilization of the iliosacral joint and the pelvic symphysis.

New sizes from August 15, 2016

#### ErgoLoc: sizes, item number and basin circumference WITHOUT pads

Regular	Large
1000 24038/0	1000 24039/0
< 99 cm	99 - 137 cm

#### ErgoLoc: sizes, item number and basin circumference WITH pads

Regular	Large
1000 24038	1000 24039

Regular	Large
< 99 cm	99-137 cm

**Aid-No. 23.11.01.0008**

**Pad accessories**

Raw material
1000 24050