

ACTIVE GENU

Knee frame orthosis



INDICATIONS

- Injuries and instability of the cruciate ligament with or without injury to anterior and posterior cruciate ligament (CI orthosis)
- Side ligament rupture or instability
- Conservative treatment of knee joint instability

MODE OF OPERATION

- Patellar subluxation "anterior and posterior drawer", stabilizes the knee joint in the sagittal area

FEATURES

- Improved fit on the inner side of the thigh
- The 4-point principle prevents pathological leverage
- Easy fitting and positioning by numbered straps
- Anatomically preformed, rigid frame geometry
- Flexion and extension stop possible
- Enclosed wedges: extension: 0 °, 10 °, 20 °, 30°, 45° | Flexion: 0 °, 10 °, 20 °, 30°, 45°, 60°, 75 °, 90 °
- Combined size as Active Genu VARIO available

Sizes, item number, thigh & calf circumference (Active Genu) *

	S	M	L	XL	XXL
Thigh	39 - 46 cm	46 - 53 cm	53 - 60 cm	60 - 67 cm	67 - 74 cm
Unter limb	28 - 34 cm	33 - 39 cm	37 - 44 cm	42 - 50 cm	48 - 57 cm
right	1000 04051	1000 04052	1000 04053	1000 04054	1000 04055
left	1000 04061	1000 04062	1000 04063	1000 04064	1000 04065

Sizes, item number, thigh & calf circumference (Active Genu VARIO) *

M / S L / M XL / L XXL / XL

Thigh	46 - 53 cm	53 - 60 cm	60 - 67 cm	67 - 74 cm
Unter limb	28 - 34 cm	33 - 39 cm	37 - 44 cm	42 - 50 cm
right	1000 04071	1000 04072	1000 04073	1000 04074
left	1000 04075	1000 04076	1000 04077	1000 04078

Aid No. 23.04.03.3052

* Circumference measured: 15 cm above or below the Centre of the patella